




- ◆ NPWT Promotes Healing in Acute and Chronic Wounds
- ◆ Optional Instillation Treatment Available for NP-800



Negative Pressure Wound Therapy System(NPWT)



NP-800/600 Lifotronic (NPWT) System

Negative Pressure Wound Therapy

NP-800/600 Features

- ◆ **Pressure Settings:**
 - From -25 mmHg to -500 mmHg
 - Providing continuous
 - Intermittent and DPC therapy
- ◆ **Canisters:** 500 ml
- ◆ **Protection with Safety alarms for:**
 - Low battery
 - Canister full
 - Full canister/line blockage
 - High/Low vacuum

Cleaning Drainage Therapy (NP-800)

- ◆ Manual cleaning (Mode D)
- ◆ Auto cleaning , Cleaning first (Mode E)
- ◆ Auto cleaning , Drainage first (Mode F)

NPWT with Instillation (NP-800)

- ◆ Reduction of bioburden
- ◆ Reduced pain at dressing change
- ◆ Moisture balance in wound bed

Back hook for easily attaching to a bed, wheelchair or IV pole

High flow rate to ensure the performance under pressure

A wide range of exudate levels can be dealt with two canister sizes

Low noise levels provide treatment in relaxation

Introducing NP-800/600

- ◆ With flexible capacity and power choices
- ◆ NP-800/600 is designed for a variety of surgical
- ◆ Chronic and traumatic wounds



Clinical Proven Benefits

- ◆ Negative Pressure Wound Therapy (NPWT) can boost wound healing via well - proved clinical mechanisms
- ◆ NPWT helps promote healing by delivering vacuum-assisted pressure blow the wound surface and forms A Closed Moist Wound Healing Environment

Macro Effects:











- ◆ Removes exudate
- ◆ Reduces edema
- ◆ Removes infectious article
- ◆ Lessens wound margins
- ◆ Increases in blood flow to the wound

Applicable Wound Types:

- ◆ Pressure ulcers
- ◆ Diabetic/neuropathic ulcers
- ◆ Venous insufficiency ulcers
- ◆ Traumatic wounds
- ◆ Post surgery and dehiscd surgical incision
- ◆ Explored fistulas
- ◆ Skin grafts

Micro Effects:

- ◆ Stimulate cell proliferation
- ◆ Reduces in bacterial load
- ◆ Increase fibroblast migration
- ◆ Promote the growth of granulation tissue

Lifotronic Product Type	Type difference	Dressing PIC	Tube PIC
NPWT-1-100A (S) NPWT-2-100A (S) NPWT-3-100A (S)	1 foam dressing(100×75×33mm) 1 sheet of film different drainage tube suite	 100A (S)	 1
NPWT-1-180A (M) NPWT-2-180A (M) NPWT-3-180A (M)	1 foam dressing(180×125×33mm) 2 sheets of film different drainage tube suite	 180A (M)	 2
NPWT-1-180B (M thin) NPWT-2-180B (M thin) NPWT-3-180B (M thin)	1 thin foam dressing(180×125×16mm) 2 sheets of film different drainage tube suite	 180B (M thin)	 2
NPWT-1-260A (L) NPWT-2-260A (L) NPWT-3-260A (L)	1 foam dressing(260×150×33mm) 2 sheets of film different drainage tube suite	 260A (L)	 3
NPWT-1-260B (L thin) NPWT-2-260B (L thin) NPWT-3-260B (L thin)	1 thin foam dressing(260×150×16mm) 2 sheets of film different drainage tube suite	 260B (L thin)	 3

NP-100

Lifotronic (NPWT) System


Negative Pressure Wound Therapy

Introducing NP-100 easy to go NPWT unit

- ◆ With stylish design and quiet operation, patients can do mobile treatment at home comfortably.
- ◆ Sleek design, tiny and light weight
- ◆ Compact size offers superb portability
- ◆ Integrated 300ml canister
- ◆ Screen lock function to avoid false operation
- ◆ Quiet treatment with Easy to use controls

NP-100 Features

- ◆ **Pressure Settings:**
Canisters: 300m
- ◆ **Safety Alarms:**
High or low vacuum
Canister full
Low battery
Negative Pressure Pump life cycle ends





NP-100 Specs:

Description	NP-100
Size (HxWxD)	60×130×160 mm
Weight	0.6 kg without the bag
Vacuum Range	-50 to -200 mmHg, -50, -75, -100, -125, -150, -175, -200
Available Modes	Continuous, Intermittent
Battery	12 Hours, 5200mAh/7.4V
Recharge Times	6 Hours
Auto Lockout	Yes, display lock out after 10 minutes
Canister	300 ml
Solidifier with canister	Yes
Display	Digital
Antibacterial Filter	Yes
Patient lock out	Yes
IV Pole attachment	Via Strap
Bed Attachment	Via Strap
Carry Bag	Yes
Alarms	High / Low Vacuum, Canister full, Low Battery, Negative Pressure Pump life cycle ends
Alarm Mute	Yes
Tubing Clamps	Yes
Power Input	Internal Lion Battery or 110-240 VAC



Lifotronic (NPWT) System

NP-800/600 Specs

Description	NP-800/600
Size (HxWxD)	220×260×200 mm
Weight	4.33 kg
Vacuum Range	-25 to -500 mmHg, -25, -50, -75, -100, -125, -150..... -500
Available Modes	NP-600: Continuous, Intermittent and DPC
	NP-800: Continuous, Intermittent and DPC; Manual cleaning; Auto cleaning , Cleaning first; Auto cleaning , Drainage first
Battery	6 Hours, 5200 mAh/7.4V
Free air flow rate	≥ 8 L/min
Recharge Times	6 Hours
Auto Lockout	Yes, display lock out after 10 minutes or Manually Lock
Canister	500 ml, marking with 50 ml increments
Solidifier with canister	Yes
Display	Digital touch screen
Antibacterial Filter	Yes
Patient lock out	Yes
IV Pole attachment	Via Strap
Bed Attachment	Via Strap
Carry Bag	Yes
Bluetooth	Yes, for mobile app monitor and control
Alarms	High / Low Vacuum, Canister full, Low Battery, Negative Pressure Pump life cycle ends
Alarm Mute	Yes
Tubing Clamps	Yes
Power Input	Internal Lion Battery or 110-240 VAC



Distributore per l'Italia



SUNMEDICAL S.r.l.
 Via Modigliani, 45 - 20054 Segrate (MI)
T +39 02 87282900
F +39 02 87282910
E info@sunmedical.it
W www.sunmedical.it

