

**pfm**medical

Quality and Experience


---

**Peripherally Inserted  
Central venous Catheter  
pfm Micro Introducer Kit**

---

The pfm Micro Introducer Kit is an accessory for the insertion of PICC lines. The high quality components ensure optimal and safe handling.

[www.pfmmedical.com](http://www.pfmmedical.com)



The **pfm Micro Introducer Kit** is an accessory for the insertion of PICC lines. The high quality components ensure optimal and safe handling.

### Benefits

#### Peelable sheath introducer

The perfect fitting of the dilator and the peelable sheath provides a smooth transition during the insertion of the sheath through the skin into the vessel.

#### Echogenic needle

The echogenic needle allows a good visualization of the needle-bevel via ultrasound to ensure proper puncture of the vessel.

#### Nitinol guidewire

The 40 cm long Nitinol guidewire is anti-kinking and ensures a simple and smooth introduction of the device.

### Technical Data

#### Set components:

- ▶ Nitinol guidewire 40 cm
- ▶ Peelable sheath introducer
- ▶ Echogenic introducer needle 21 G / 7 cm

### Ordering Information

REF	Size Peelable sheath	PU
P090700011	3.0 F (1.0 mm)	5
P090700012	4.0 F (1.3 mm)	5
P090700013	5.0 F (1.7 mm)	5
P090700014	6.0 F (2.0 mm)	5
P090700015	7.0 F (2.3 mm)	5

### Contact

Should you have any questions our Customer Solutions Team  is happy to advise you.

service@pfmmmedical.com

T +49 (0)2236 9641-220

F +49 (0)2236 9641-51

[www.pfmmmedical.com](http://www.pfmmmedical.com) - Your information source for the products of pfm medical ag

**pfm**medical

pfm medical ag  
Wankelstraße 60  
50996 Köln, Germany

Certified by  
DIN EN ISO 13485

Manufacturer:  
pfm medical usa Inc.  
1815 Aston Ave, Suite 106  
Carlsbad, CA 92008, USA